

ADORE

Supporting Health



CATALOG 2021

Pain Management Aids

Rehabilitation Aids

Immobilizers

Traction Kits



ADORE

Supporting Health

Our mission is to make high quality orthopaedic support braces, fracture aids and pain management aids easily accessible because we believe that everyone deserves good health.

We strictly follow international standards and source only the highest quality materials to develop our products.

Our commitment to supporting good health for everyone drives every facet of our business and partnerships. Innovation, dedication and strong relationships with our distribution partners are the core values that define Adore Medical Private Limited.



CONTENTS



06

ABDOMEN & SPINE SUPPORT

Abdominal Belt
Tummy Trimmer
Maternity Belt
Ash Brace
Rib Belt
Sternal Support
Lumbo Sacral Support
Contoured L.S. Support
Dorso Lumbar Brace
Hernia Belt
Waist Belt (Neoprene)
Scrotal Support
Posture Corrector
Back Rest
Coccyx Cushion



28

ARM & ELBOW SUPPORT

Pouch Arm Sling (Tropical)
Pouch Arm Sling (Baggy)
Elbow Support
Tennis Elbow Support
Forearm Splint
Tourniquet



16

CERVICAL SUPPORT

Cervical Collar Hard
Cervical Collar Soft
Cervical Collar Soft with Support
Cervical Pillow (Memory Foam)
Cervical Pillow (PU FOAM)
Cervical Traction Kit (With Weight Bag)
Philadelphia Collar (Cervical Orthosis)



22

SHOULDER & CLAVICLE SUPPORT

Shoulder Support (Neoprene)
Elastic Shoulder Immobilizer
Universal shoulder Immobilizer
Clavicle Brace
Clavicle Brace With Buckle



34

WRIST FINGERS & THUMB SUPPORT

Wrist Support (Neoprene)
Wrist Support Double Lock
Wrist Support with Thumb
Wrist Support With Thumb (Neoprene)
Wrist & Forearm Splint
Elastic Wrist Splint
Finger Cot Splint
Finger Extension Splint
Frog Splint
Mallet Finger Splint
Thumb Spica Splint
Exercise Gel Ball (Soft/Medium/Hard)
Weight Cuff



42

PELVIC & HIP SUPPORT

Pelvic Binder
Pelvic Traction Kit
Leg Traction Brace



58

FOOT, HEEL & TOES SUPPORT

Foot Drop Splint
Arch Support
Bunion Splint
Toe Separator
Insole
Heel Cushion (silicone)
Metatarsal Pad
Cast Shoe



64

HOT & COLD THERAPY AIDS

Cool Eye Pack
Hot & Cold Pack



46

THIGH, KNEE, ANKLE & CALF SUPPORT

Skin Traction (PUF Liner)
Thigh Support
Ankle Brace
Ankle Binder
Anklet
Knee Cap
Knee Cap with Open Patella
Knee Cap with Rigid Hinge
Knee Wrap (Neoprene)
Knee Wrap Hinged (Neoprene)
Elastic Knee Support
Patellar Support
Compression Stocking Below Knee
Compression Stocking Mid Thigh
Knee Immobilizer (Short type)
Knee Immobilizer (Long type)
Knee Immobilizer (22")
R.O.M. Knee Brace
Calf Support



ADORE
Supporting Health



ABDOMEN & SPINE SUPPORT

ABDOMINAL BELT

TUMMY TRIMMER

MATERNITY BELT

ASH BRACE

RIB BELT

STERNAL SUPPORT

LUMBO SACRAL SUPPORT

CONTOURED L.S. SUPPORT

DORSO LUMBAR BRACE

HERNIA BELT

WAIST BELT (NEOPRENE)

SCROTAL SUPPORT

POSTURE CORRECTOR

BACK REST

COCCYX CUSHION



AD-101




ABDOMINAL BELT

Our abdominal belt is perfect for postoperative conditions of the abdomen and post-delivery umbilical and ventral hernias. It features easy-to-use hook & loop closures that ensures a proper fit. High quality elastic materials with a breathable and moisture-wicking weave ensure a comfortable fit.

INDICATIONS:

Postoperative conditions of the abdomen.
Post-delivery umbilical and ventral hernia.

| SIZE | INCHES | CM |  |
|-------|--------|---------|--|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |
| XXXXL | 52-56 | 130-140 | |



AD-107




TUMMY TRIMMER

Optimal compression of our Tummy Trimmer binds surgical incisions and helps heal abdomen and chest. Ideal to relieve pain associated with movement or cough and support during workout. Porous three-panel webbing and a broad hook and loop tape panel provide optimal compression to tone up abdominal muscles post surgery or pregnancy.

INDICATIONS:

| | |
|----------------------------------|------------------------------------|
| Postural deformity. | Mild thoraces and lumbar injuries. |
| Osteoporosis of the spine. | Degenerative Disc Disease |
| Postoperative care of the spine. | Herniated Disc or PIVD. |

| SIZE | INCHES | CM |  |
|-------|--------|---------|---|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |
| XXXXL | 52-56 | 130-140 | |

AD-111




MATERNITY BELT

This maternity belt provides relief of hip and pelvic pain through sufficient stretching and gentle compression to reduce bladder pressure and spread the baby's weight. It supports weak abdominal muscles during pregnancy as well as after pregnancy.

INDICATIONS:

Postoperative conditions of the abdomen.
Post-delivery umbilical and ventral hernia.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 34-38 | 86-97 | |
| M | 38-42 | 97-107 | |
| L | 42-46 | 107-117 | |
| XL | 46-50 | 117-127 | |
| XXL | 50-54 | 127-137 | |



AD-108



ASH BRACE

Our Ash Brace is aerodynamically shaped with telescoping, aluminum section, and malleable horizontal main body to form a cross bar. Knob and screw mechanism allows two-way size adjustment of the brace on the vertical side. Offers pressure relief in thoracic and lumbar spine and moderate immobilization.


INDICATIONS:

Osteoporotic Bone Collapse.

Ankylosing Spondylitis.

Conservative management of Vertebral Bone Tuberculosis.

Postoperative, post-traumatic, or chronic spinal pathology. Complications of thoracolumbar junction which needs moderate Immobilization.

| SIZE | FEET | CM |  |
|-------|---------|---------|---|
| Short | 5-6 | 152-183 | |
| Long | Above 6 | > 183 | |



AD-104




RIB BELT

Designed with a splinting pad for comfort, our Rib Belt is useful in aiding the healing process in rib fracture, relieving chest pain, and intercostal muscle strain. Features an easy hook & loop closure and a breathable, moisture-wicking weave for comfort and quick application and removal.

INDICATIONS:

- Postoperative conditions of the abdomen.
- Post-delivery.
- Umbilical and ventral hernia.

| SIZE | INCHES | CM |  |
|------|--------|---------|--|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |



AD-105




STERNAL SUPPORT

Designed for comfort with an easy hook & loop closure, moisture-wicking fabric and splinting pad. It is also lightweight to allow for maximum breathability and prevent sweating discomfort. Helps keep the rib cage intact to prevent any additional risk of injury and easy wearing with an adjustable hook & loop closure.

INDICATIONS:

- Post cardiothoracic surgeries.
- Fracture of the sternum.
- Pain due to chest injury.
- Rib fractures.
- Osteoporosis.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |

AD-102




LUMBO SACRAL SUPPORT

Our lumbo-sacral support is used to comfortably maintain correct posture and support the lower back while sitting, walking or driving. It is ideal for use during postoperative care of the lumbar spine, lumbar disc injury, degenerative disc disease, herniated disc/PVD, and for low back pain, strain and sprains.

INDICATIONS:

Postoperative conditions of the abdomen.
Post-delivery umbilical and ventral hernia.

| SIZE | INCHES | CM |  |
|-------|--------|---------|--|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |
| XXXXL | 52-56 | 130-140 | |



AD-103




CONTOURED L.S. SUPPORT

This Contoured Lumbo-Sacral Support is ideal for postoperative care of lumbar spine and for lower back pain, strain or sprains. Breathable, moisture-wicking fabric makes it comfortable for extended wear while sitting, walking or driving. Hook & loop closure and side pulls ensure proper fit, added support, and quick and easy application.

INDICATIONS:

Postoperative care of the lumbar spine.
Low back pain.
Low back sprains/strains.
Lumbar disc injury.

Degenerative Disc Disease.
Herniated Disc or PIVD.

| SIZE | INCHES | CM |  |
|-------|--------|---------|---|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |
| XXXL | 48-52 | 120-130 | |
| XXXXL | 52-56 | 130-140 | |



AD-106




DORSO LUMBAR BRACE


Features anatomically contoured splints for precise immobilization and clavicle straps with cotton stockinet to correct shoulder posture alignment. Thoracic straps provide additional support. Easy hook and loop closures on abdomen elastic panel for proper fit and easy application and removal.

INDICATIONS:

Postural deformity.
Osteoporosis of the spine.
postoperative spine care.
Mild lumbar and thoraces injury.

Degenerative Disc Disease.
Herniated Disc or PIVD.

| SIZE | FEET | CM |  |
|-------|---------|---------|---|
| Short | 4.6-5.6 | 140-171 | |
| Long | 5.6-6.6 | 171-201 | |

| SIZE | INCHES | CM |  |
|-----------|--------|---------|--|
| Universal | 28-44 | 70-110 | |
| Special | 44-52 | 110-130 | |

AD-109




HERNIA BELT

Designed for a comfortable and supportive fit our hernia belt is lightweight and features adjustable pelvic and leg straps. Soft, dermophilic material is gentle on sensitive skin and deters allergies and rashes.

INDICATIONS:

Reducible inguinal or scrotal hernia.
Post hernia repair surgery.
Post-delivery umbilical and ventral hernia.
Post abdominal hernia repair surgery.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |

AD-112




WAIST BELT (NEOPRENE)

Our wrap-around design waist belt is made with neoprene bonded fabric for four-way stretchability and maximum support. Material retains body heat to provide quick and consistent pain relief. Provides abdominal muscle support during walking, running, aerobics, and other sports activities.

INDICATIONS:

Weak abdominal muscles.
Swelling & stiffness in abdominal muscles.
Postoperative, post-delivery support.

| SIZE | INCHES | CM |  |
|-----------|--------|--------|--|
| Universal | 28-50 | 70-125 | |



AD-1103




SCROTAL SUPPORT

Fabricated extra-soft and hypoallergenic materials, our scrotal support is well-ventilated to allow for ease of use and extended wear comfort. Features a separated penile compartment and adjustable waist band.

INDICATIONS:

Palliative care.
Treatment for postoperative care.
Tenderness, pain or swelling in scrotal area.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |



AD-110



POSTURE CORRECTOR

Designed to correct slouching and improve posture. Features breathable fabric with cotton lining for comfortable and easy extended wear. Can be worn while working, walking and exercising.

INDICATIONS:

- Back pain.
- Lumbar pain.
- Spine misalignment.
- Poor posture.

| SIZE | INCHES | CM | |
|------|--------|---------|--|
| XS | 22-28 | 56 - 71 | |
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |



AD-1101



BACK REST

Promotes good posture and correct spinal alignment by allowing the spine to follow natural curve. Alleviates back pain and provides support for lower back.

INDICATIONS:

- Postoperative discomfort.
- Post-injury pain.
- Poor posture.
- Lower back pain.
- Pregnancy discomfort.

| SIZE | CM | |
|-----------|---------|--|
| Universal | 152-192 | |

AD-1102



COCCYX CUSHION

The high-density foam of our coccyx cushion features a U-shaped cut that gently suspends your coccyx to relieve pressure on your spinal cord. Contoured to improve posture and help sustain a cervical curve that is optimal for improved blood flow.

INDICATIONS:

- Pudendal nerve entrapment pain.
Chronic pelvic pain syndrome.
Rectal surgery.
Prostatitis.
Anal surgery.
- Coccyx fracture.
Hemorrhoids.
Interstitial Cystitis.
Pelvic floor pain.
Pressure sores.

| SIZE | INCHES (LxBxH) | CM (LxBxH) |
|-----------|-------------------|---------------|
| Universal | 16.4x16.5x2.7 | 41.7x42x7 |



ADORE
Supporting Health



CERVICAL SUPPORT AIDS

CERVICAL COLLAR HARD

CERVICAL COLLAR SOFT

CERVICAL COLLAR SOFT WITH SUPPORT

CERVICAL PILLOW (MEMORY FOAM)

CERVICAL PILLOW (PU FOAM)

CERVICAL TRACTION KIT

PHILADELPHIA COLLAR



AD-205




CERVICAL COLLAR HARD

Our hard Cervical Collar is made of high-quality polyethylene with padded edges moulded with Ethafoam. Ethafoam padding offers high cushioning for improved comfort and wear. Air ventilation holes allow proper air circulation and comfort. Collar height can be adjusted accurately to fit any neck height.

INDICATIONS:

Cervical Spondylitis.
Cervical Disc Disease.
Fracture of cervical spine.
Whiplash injuries.

| SIZE | INCHES | CM |  |
|------|---------|----------|--|
| S | 3.0-3.5 | 7.5-8.75 | |
| M | 3.5-4.0 | 8.75-10 | |
| L | 4.0-4.5 | 10-11.25 | |



AD-202




CERVICAL COLLAR SOFT

Made of soft density foam for gentle and cushioned neck support. Our soft Cervical Collar is fabricated with a soft, smooth cotton stockinet and has hook & loop closers to allow proper fit. This is a durable cervical collar for long-term use.

INDICATIONS:

| | |
|--------------------------|----------------------------|
| Cervical Trauma. | Strain of Neck Muscles. |
| Neck Pain and Stiffness. | Whiplash Recovery. |
| Excessive Mobility. | Degenerative Disc Disease. |
| Rheumatoid Arthritis. | |

| SIZE | INCHES | CM |  |
|------|---------|-------------|---|
| S | 3.0-3.5 | 7.5-8.75 | |
| M | 3.5-4.0 | 8.75-10 | |
| L | 4.0-4.5 | 10-11.25 | |
| XL | 4.0-5.0 | 11.25-12.50 | |

AD-201




CERVICAL COLLAR SOFT WITH SUPPORT

Made of firm density foam covered with a soft, smooth cotton stockinet splint inside for better neck support. Eyelets provide air ventilation and hook & loop closers allow proper fit. This is a durable cervical collar for long-term use.

INDICATIONS:

| | |
|--------------------------|----------------------------|
| Cervical Trauma. | Strain of neck muscles. |
| Neck Pain and stiffness. | Whiplash recovery. |
| Excessive mobility. | Degenerative Disc Disease. |
| Rheumatoid Arthritis. | |

| SIZE | INCHES | CM |  |
|------|---------|-------------|--|
| S | 3.0-3.5 | 7.5-8.75 | |
| M | 3.5-4.0 | 8.75-10 | |
| L | 4.0-4.5 | 10-11.25 | |
| XL | 4.5-5.0 | 11.25-12.50 | |



AD-203




CERVICAL PILLOW (MEMORY FOAM)

Made of memory foam to provide proper support to the neck region while relaxing or sleeping. This Cervical Pillow is covered with a soft, smooth, and skin-friendly cotton fabric.

INDICATIONS:

| | |
|--------------------------------------|--------------------|
| Temporomandibular disorders. | Whiplash injuries. |
| Cervical trauma. | Arthritis. |
| Morning hand stiffness and swelling. | Neck pain. |

| SIZE | INCHES | CM |  |
|-----------|--------|----------|---|
| Universal | 3-5 | 7.5-12.5 | |



AD-206




CERVICAL PILLOW (PU FOAM)

Provides sustained support to cervical vertebrae during sleep. This cervical pillow is designed to correct and retain cervical alignment.

INDICATIONS:

- Temporomandibular disorders.
- Cervical trauma.
- Morning hand stiffness and swelling.
- Whiplash injuries.
- Arthritis.
- Neck pain.

| SIZE | INCHES | CM |  |
|-----------|--------|----------|--|
| Universal | 3-5 | 7.5-12.5 | |

AD-802



CERVICAL TRACTION KIT + WEIGHT BAG

With a width adaptation from 25mm to 50mm, our Cervical Traction Kit is developed to fit on the top of most doors. Two holding pads connected with cord loops hold the chin and the occipital area of the skull. The weight bag holds the loops of cord coming out of two sides of the head halter. Smooth surface design minimizes traction loss due to friction.

INDICATIONS:

- Cervical Spondylitis.
- Cervical Disc Disease.
- Fractures of cervical spine.
- Mechanical stress due to poor posture.

SIZE: Universal

AD-204




PHILADELPHIA COLLAR (CERVICAL ORTHOSIS)

Made of Ethafoam material to padded support our Philadelphia Collar has a trachea opening in front. The hook & loop closures enable proper fit and two-piece design provides rigid immobilization of the neck. Chin support ensures extra support to the neck.

INDICATIONS:

- Cervical trauma.
- Cervical sprain.
- Cervical Spondylitis.

| SIZE | INCHES | CM |  |
|------|---------|----------|--|
| S | 3.0-3.5 | 7.5-8.75 | |
| M | 3.5-4.0 | 8.75-10 | |
| L | 4.5-5.0 | 10-11.25 | |



ADORE
Supporting Health



SHOULDER & CLAVICLE SUPPORT

SHOULDER SUPPORT (NEOPRENE)

ELASTIC SHOULDER IMMOBILIZER

UNIVERSAL SHOULDER IMMOBILIZER

CLAVICLE BRACE

CLAVICLE BRACE
WITH BUCKLE



AD-305



SHOULDER SUPPORT (NEOPRENE)

Made from Neoprene bonded nylon fabric which is four-way stretchable with an anatomical shape to ensure perfect fitting and optimal immobilization. Retains body warmth for pain relief and faster recovery. Easily worn and can be used unobtrusively under clothing.

INDICATIONS:

- Shoulder Injury.
Rotator cuff injury.
Shoulder dislocation
Shoulder inflammation.
- Arthritic pain of the shoulder joint.
Bursitis and Tendonitis.

| SIZE | INCHES | CM | |
|-----------|--------|------------|--|
| Universal | 26-36 | 66-91.4 | |
| Special | 36-48 | 91.4-121.4 | |



AD-302



ELASTIC SHOULDER IMMOBILIZER

Our Elastic Shoulder Immobilizer is made from breathable and moisture transporting weaving fabric. Featuring a shoulder cup to reinforce immobilization, a chest band for controlling compression and application, along with humerus and forearm cuffs for extra support and balance.

INDICATIONS:

- Shoulder dislocation.
Post-surgical immobilization.
Shoulder joint Pain.
Forearm pain.
Shoulder sprain.

| SIZE | INCHES | CM | |
|------|--------|---------|--|
| S | 28-32 | 70-80 | |
| M | 32-36 | 80-90 | |
| L | 36-40 | 90-100 | |
| XL | 40-44 | 100-110 | |
| XXL | 44-48 | 110-120 | |

AD-303




UNIVERSAL SHOULDER IMMOBILIZER

Our Universal Shoulder Immobilizer comes in a flexible size to support shoulder immobilization. Featuring a baggy sling with an extra wrap-around strap and hook & loop closures for proper fit. Made from breathable and moisture transporting weaving fabric.

INDICATIONS:

Shoulder dislocation.
Post-surgical immobilization.
Shoulder joint Pain.
Forearm pain.
Shoulder sprain.

| SIZE | INCHES | CM |  |
|-----------|--------|------------|--|
| Universal | 26-36 | 66-91.4 | |
| Special | 36-48 | 91.4-121.4 | |



AD-304




CLAVICLE BRACE

Eight-shaped band design of Clavicle Brace is cushioned from soft polyurethane foam padding. It offers an elastic, ready to use application, while allowing full arm movement and immobilization of the clavicle bone.

INDICATIONS:

Clavicle fracture.
Clavicle strain.
Clavicle trauma.
Poor posture problems.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 28-34 | 70-85 | |
| M | 34-40 | 85-100 | |
| L | 40-46 | 100-115 | |
| XL | 46-52 | 115-130 | |



AD-307




CLAVICLE BRACE WITH BUCKLE

Eight-shaped band design of Clavicle Brace is cushioned from soft polyurethane foam padding. It offers adjustable straps for custom fit and compression, while allowing full arm movement and immobilization of the clavicle bone.

INDICATIONS:

- Clavicle fracture.
- Clavicle strain.
- Clavicle trauma.
- Poor posture problems.

| SIZE | INCHES | CM |  |
|------|--------|---------|--|
| S | 28-34 | 70-85 | |
| M | 34-40 | 85-100 | |
| L | 40-46 | 100-115 | |
| XL | 46-52 | 115-130 | |



ADORE

Supporting Health



ADORE
Supporting Health



ARM & ELBOW SUPPORT

POUCH ARM SLING (TROPICAL)

POUCH ARM SLING (BAGGY)

ELBOW SUPPORT

TENNIS ELBOW SUPPORT

FOREARM SPLINT

TOURNIQUET



AD-301




POUCH ARM SLING (TROPICAL)

Our Pouch Arm Sling provides comfortable support to the wrist and arm with an adjustable strap with a buckle for proper immobilization. It offers a breathable, moisture-wicking weave for comfort and quick application.

INDICATIONS:

- Arm fracture.
- Wrist fracture.
- Arm & Wrist sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast.

| SIZE | INCHES | CM |  |
|------|-----------|---------|--|
| XS | 19.7-23.6 | 50-60 | |
| S | 27.6-31.5 | 70-80 | |
| M | 31.5-35.4 | 80-90 | |
| L | 34.4-39.4 | 90-100 | |
| XL | 39.4-43.3 | 100-110 | |



AD-306




POUCH ARM SLING (BAGGY)

Our baggy Pouch Arm Sling provides light-weight support to the wrist and arm through a thumb cradle, and an adjustable strap with a buckle for proper immobilization. It offers a breathable, moisture-wicking weave for comfort and quick application.

INDICATIONS:

- Arm Fracture.
- Wrist Fracture.
- Arm & Wrist Sprains.
- Arm/Shoulder immobilization.
- Post-surgery/Post-cast rehabilitation.

| SIZE | INCHES | CM |  |
|------|-----------|---------|---|
| XS | 19.7-23.6 | 50-60 | |
| S | 27.6-31.5 | 70-80 | |
| M | 31.5-35.4 | 80-90 | |
| L | 34.4-39.4 | 90-100 | |
| XL | 39.4-43.3 | 100-110 | |

AD-504




ELBOW SUPPORT

Our easy-to-wear compression support for the elbow is made from sturdy, fully stretchable fabric. Comfortable and moisture-transporting material makes it useful to wear in any weather.

INDICATIONS:

- Forearm/Elbow pain.
Post Elbow
Arthroscopy.
Sports and workout
stress.
- Post-cast
rehabilitation.
Postoperative
weakness.
Arthritis / Bursitis.

| SIZE | INCHES | CM |  |
|------|--------|-------|---|
| S | 7-12 | 18-25 | |
| M | 10-12 | 25-30 | |
| L | 12-14 | 30-35 | |
| XL | 14-16 | 35-40 | |



AD-503




TENNIS ELBOW SUPPORT

Our Tennis Elbow support has adjustable straps and a reverse fastener for a strong grip and custom compression. Comfortable and moisture-transporting material makes it useful to wear anytime.

INDICATIONS:

- Tennis/Golf Elbow.
Tendonitis.
Swelling.

| SIZE | INCHES | CM |  |
|------|--------|-------|---|
| S | 7-12 | 18-25 | |
| M | 10-12 | 25-30 | |
| L | 12-14 | 30-35 | |
| XL | 14-16 | 35-40 | |



AD-507




FOREARM SPLINT

Our Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb. The triple-layered bonded material with the foam padding midlayer provides soft cushioning and makes it comfortable to wear.

INDICATIONS:

Post-cast/Post-surgery protection and rehabilitation.
Conservative management of Kienbock's disease.
Rheumatic disease of the wrist.
Arthritic hand positioning.

| SIZE | INCHES | CM |  |
|-----------|--------|-------|--|
| Universal | 7.2-12 | 18-30 | |



AD-1201



TOURNIQUET

(UPPER & LOWER EXTREMITIES)

Our Tourniquet is padded with soft foam on the inner surface and constructed with a lasting hook & loop strap with a tough polypropylene fastener.

INDICATIONS:

Stop constant bleeding.
Restrain the circulation of blood to an extremity for a period of time.

| SIZE |
|-----------|
| Universal |

ADORE

Supporting Health



ADORE
Supporting Health



WRIST FINGERS & THUMB SUPPORT

WRIST SUPPORT (NEOPRENE)

WRIST SUPPORT DOUBLE LOCK

WRIST SUPPORT WITH THUMB

WRIST SUPPORT WITH THUMB (NEOPRENE)

WRIST & FOREARM SPLINT

ELASTIC WRIST SPLINT

FINGER COT SPLINT

FINGER EXTENSION SPLINT

FROG SPLINT

MALLET FINGER SPLINT

THUMB SPICA SPLINT

EXERCISE GEL BALL

WEIGHT CUFF



AD-508




WRIST SUPPORT (NEOPRENE)

Created out of triple-layered neoprene. The top layer features continuous loops with rugged nylon fabric in the inner layer. Anatomical thumb slot designed for limited abduction of the thumb.

INDICATIONS:

- Tendonitis.
- Tenosynovitis.
- Arthritic care of wrist.
- Sprains and strains.
- Dislocation of the wrist joint.
- Postoperative and post-cast care.
- Prophylaxis for exercise and sports.

| SIZE | INCHES | CM |  |
|-----------|---------|-----------|--|
| Universal | 4.5-9.5 | 11.4-23.8 | |

AD-501




WRIST SUPPORT DOUBLE LOCK

Designed with adjustable straps, moisture-wicking fabric and durable nylon reinforcement. Perfect companion for sport activities and for comfort from swelling and sprains of the wrist.

INDICATIONS:

- Sprains and strains.
- Mild Osteoarthritis.
- Swelling & inflammation.
- Prophylaxis for exercise and sports.

| SIZE | INCHES | CM |  |
|------|---------|-----------|---|
| S | 4.9-6.1 | 12.5-15.5 | |
| M | 6.1-7.3 | 15.5-18.5 | |
| L | 7.3-8.5 | 18.5-21.5 | |

AD-502




WRIST SUPPORT WITH THUMB

Designed with adjustable wrap-around straps, moisture-wicking fabric and durable nylon reinforcement. Supports upper wrist area and has a thumb opening for strong grip and compression.

INDICATIONS:

Sprains and strains .
Mild Osteoarthritis.
Swelling & inflammation.
Tendonitis.

| SIZE | INCHES | CM |  |
|-----------|---------|-----------|--|
| Universal | 4.5-9.5 | 11.4-23.8 | |



AD-509




WRIST SUPPORT WITH THUMB (NEOPRENE)

Neoprene based design with adjustable wrap-around straps, moisture-wicking fabric and durable nylon reinforcement. It retains body warmth to relieve pain and accelerate recovery.

INDICATIONS:

Ligamentous laxity of the thumb.
Median nerve injury.
Skier's thumb.
Carpometacarpal joint injury.
Post operative and post-cast care.
Prophylaxis for exercise and sports.

| SIZE | INCHES | CM |  |
|-----------|---------|-----------|---|
| Universal | 4.5-9.5 | 11.4-23.8 | |



AD-505




WRIST & FOREARM SPLINT

Our Wrist & Forearm Splint is made of lasting nylon sleeve supported from within with pliable aluminum, restricts wrist movement, but allows the full extent of motion to the fingers and thumb. The triple-layered bonded material with the foam padding midlayer provides soft cushioning and makes it comfortable to wear.

INDICATIONS:

- Post-cast management.
- Post-surgical protection and rehabilitation.
- Conservative management of Kienbock's disease.
- Rheumatic disease of the wrist.
- Arthritic hand positioning.

| SIZE | INCHES | CM |  |
|------|------------|-----------|--|
| S | 7.2-8.5 | 18.3-21.6 | |
| M | 8.5-9.8 | 21.6-24.9 | |
| L | 9.8 - 11.1 | 24.9-28.2 | |
| XL | 11.1-12.4 | 28.2-31.5 | |



AD-506




ELASTIC WRIST SPLINT

Our athletic-inspired elastic splint is a lightweight wrist support with double straps and anatomical design for rigid immobilization, while maintaining the wrist in a functional position.

INDICATIONS:

- Carpal Tunnel Syndrome.
- Post-surgery.
- Post-cast management.
- Wrist sprains and strains.
- Tendonitis.

| SIZE | INCHES | CM |  |
|------|------------|-----------|---|
| S | 7.2-8.5 | 18.3-21.6 | |
| M | 8.5-9.8 | 21.6-24.9 | |
| L | 9.8 - 11.1 | 24.9-28.2 | |
| XL | 11.1-12.4 | 28.2-31.5 | |

AD-601




FINGER COT SPLINT

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort. Loop tape adjuster provides fastening to apply in a fitting position.

INDICATIONS:

- Finger injury.
- Post finger Fracture.
- Postoperative care.
- Sprains and bruises of fingers.

| SIZE | INCHES | CM |  |
|------|--------|--------|--|
| S | < 2 | < 5 | |
| M | 2-3 | 5 -7.5 | |
| L | 3-4 | 7.5-10 | |



AD-604




FINGER EXTENSION SPLINT

Made of aluminum exterior frame for protective and perfect immobilization and soft foam inner liner for comfort. The built-in velcro tape on the finger and the wrist ensures a tight fit, excellent grip and easy application.

INDICATIONS:

- Metacarpal injury.
- Fracture.
- Post-surgery.
- Sprains and bruises of fingers.

| SIZE | INCHES | CM |  |
|------|--------|--------|---|
| S | < 2 | < 5 | |
| M | 2-3 | 5 -7.5 | |
| L | 3-4 | 7.5-10 | |



AD-602




FROG SPLINT

Made of light-weight but rigid aluminum frame to rehabilitate and align the post ligamentary damage of the Interphalangeal joints, our Frog Splint has soft foam inner liner for comfort and aluminum wings to adjust for a perfect fit.

INDICATIONS:

- Interphalangeal (IP) joint Injury.
- Finger Injury.
- Finger sprain.
- Postoperative care.

| SIZE | INCHES | CM |  |
|------|--------|--------|--|
| S | < 2 | < 5 | |
| M | 2-3 | 5 -7.5 | |
| L | 3-4 | 7.5-10 | |

AD-605




MALLET FINGER SPLINT

Ethafoam liner for cushioning and also suitable for post-operative care. Malleable aluminum frame can be adjusted alignment and immobilization. Can be easily worn for most sports and daily activities.

INDICATIONS:

- Bruising of the dorsum of the DIPJ.
- Mallet finger deformity.
- May help manage the extensor tendon injury.
- Minor fractures of the distal phalanges.
- Distal phalange post-fracture care.
- Distal interphalangeal (DIP) joint conditions.
- Swollen fingertip.
- Post-surgical care.

| SIZE | INCHES | CM |  |
|-----------|--------|--------|---|
| Universal | < 2-4 | < 5-10 | |

AD-603




THUMB SPICA SPLINT

A flexible and comfortable splint made of malleable material for easy application and removal. It allows free movement of the hand and fingers while supporting the thumb. Its skin-friendly breathable fabric allows moisture transport and comfort.

INDICATIONS:

| | |
|----------------------------------|-----------------------------|
| Bursitis. | Thumb pain symptoms. |
| Sprains of thumb. | CMC & MCP joint conditions. |
| Avulsion fractures of the thumb. | |
| Soft tissue injuries. | |

| SIZE | INCHES | CM |  |
|-----------|--------|-----------|--|
| Universal | 5-9.5 | 12.5-23.8 | |



AD-902



EXERCISE GEL BALL

● **SOFT** ● **MEDIUM** ● **HARD**

Anti-stress exercise gel ball in three colours depending on hardness. Good for building strength of fingers, palm and wrist.

INDICATIONS:

| | |
|----------------|---------------------------------|
| Post fracture. | Weak forearm and grip strength. |
| Stress. | |

AD-901



WEIGHT CUFF

Anti-stress exercise gel ball in three colours depending on hardness. Good for building strength of fingers, palm and wrist.

INDICATIONS:

| | |
|----------------|---------------------------------|
| Post fracture. | Weak forearm and grip strength. |
| Stress. | |

| SIZE | ½ Kg | 1Kg | 2Kg |
|--|------|-----|-----|
| Variations of Weight Cuff fit all wrists and ankles. | | | |



ADORE
Supporting Health



PELVIC & HIP SUPPORT

PELVIC BINDER

PELVIC TRACTION KIT (WITH WEIGHT BAG)

LEG TRACTION BRACE

PRODUCT CATALOG 2021



AD-801




PELVIC BINDER

Our Pelvic Binder is made with ruggedness for extended and repetitive uses. With a large hook loop area, it can accommodate customized compression and adaptable sizing. Inclined traction straps ensure the directional force onto the vertebrae in the lower part of the spine. This can be easily worn while under care at home or in hospitals.

INDICATIONS:

- Lumbar spondylitis.
- Postural deformities.
- Lower back pain.
- Vertebral fractures.
- Postural fatigue.
- Mechanical stress due to Poor posture.

| SIZE | INCHES | CM |  |
|------|--------|---------|--|
| S | 32-36 | 80-90 | |
| M | 36-40 | 90-100 | |
| L | 40-44 | 100-110 | |
| XL | 44-48 | 110-120 | |
| XXL | 44-52 | 120-130 | |

AD-803




PELVIC TRACTION KIT (WITH WEIGHT BAG)

A versatile traction kit that can be used in all kinds of tractions and on all kinds of beds. The traction strap braces the pelvic area and pulls the vertebral column. **Kit includes** - Traction Pulley Bracket, Pelvic Binder, Spreader Bar, Traction Cord & Weight Bag.

INDICATIONS:

- Lumbar Spondylitis.
- Postural deformities.
- Lower back pain.
- Vertebral fractures.
- Postural fatigue.
- Mechanical stress due to Poor posture.

| SIZE | INCHES | CM |  |
|------|--------|---------|---|
| S | 32-36 | 80-90 | |
| M | 36-40 | 90-100 | |
| L | 40-44 | 100-110 | |
| XL | 44-48 | 110-120 | |

AD-804



LEG TRACTION BRACE

Four sets of nylon based buckle closure with an additional locking mechanism to ensure a secure fit. Strong straps with metallic hoops enable easy binding to the traction system. Made from textured foam for comfort and secure grip.

INDICATIONS:

- Traction to relieve pain in the lower back when pelvic traction is not possible or not desired.
- Hip/Sacral pain.
- Shaft fractures of the femur.

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 9.8-11.3 | 25-29 | |
| M | 11.3-12.8 | 29-33 | |
| L | 12.8-14.3 | 33-36 | |
| XL | 14.3-15.8 | 36-40 | |



ADORE
Supporting Health



THIGH, KNEE, ANKLE & CALF SUPPORT

SKIN TRACTION (PUF LINER)

THIGH SUPPORT

ANKLE BRACE

ANKLE BINDER

ANKLET

KNEE CAP

KNEE CAP WITH OPEN PATELLA

KNEE CAP WITH RIGID HINGE

KNEE WRAP (NEOPRENE)

KNEE WRAP HINGED (NEOPRENE)

ELASTIC KNEE SUPPORT

PATELLAR SUPPORT

COMPRESSION STOCKING BELOW KNEE

COMPRESSION STOCKING MID THIGH

KNEE IMMOBILIZER (SHORT TYPE)

KNEE IMMOBILIZER (LONG TYPE)

KNEE IMMOBILIZER (22")

R.O.M. KNEE BRACE

CALF SUPPORT



AD-805



SKIN TRACTION (PUF LINER)

Our Skin Traction kit is made with foam liner designed to ensure there is no friction damage to the skin and no slippage during traction. It includes a durable crepe-bandage to secure the foam liner. It is easy to apply with a built-in cord and can be used for the Upper & Lower Extremities.

INDICATIONS:

Need for indirect traction to extremities.
Pains & spasms.

| SIZE | INCHES | CM |
|-------|--------|------|
| Child | < 12 | < 30 |
| Adult | > 12 | > 30 |

AD-416



THIGH SUPPORT

The weaving method for our Thigh Support provides seamless stretching stockings with a tapered top to bottom design to hold in a natural shape. The high-quality nylon and spandex yarn used for compression is woven on state-of-the-art circular looms. The flexible thread runs spirally in the pattern of nylon locks to ensure free movement and create uniform compression even on uneven areas around the limbs.

INDICATIONS:

Hamstring strain. Arthritis.
Muscle spasms or Stiffness or trauma.
cramping. Pain & inflammation.
Sports Injury.

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |

AD-412




ANKLE BRACE

Easy to wear and remove brace made of rigid polypropylene material for strong immobilization. Further reinforcement of proper fit and comfort via a lace pull mechanism. Skin-friendly material of the brace ensures longer duration application and proper ventilation to avoid skin irritation.

INDICATIONS:

Mild ankle sprains
Recurring ankle injury.
Tarsal tunnel syndrome.
Swollen ankle.
(Not indicated for severe ankle instabilities)

| SIZE | INCHES | CM |  |
|------|----------|-------|--|
| S | 7.2-8.4 | 18-21 | |
| M | 8.4-9.6 | 21-24 | |
| L | 9.6-10.8 | 24-27 | |
| XL | 10.8-12 | 27-30 | |



AD-401




ANKLE BINDER

The figure of eight design makes it comfortably wearable, even with shoes. Spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing. High denier count provides durability and sturdiness for extended use.

INDICATIONS:

Supination injuries. Chronic ligament
Post-supination instability in the upper
injuries. & lower ankle.
Syndesmotic ligament
injury.

| SIZE | INCHES | CM |  |
|------|----------|-------|---|
| S | 7.2-8.4 | 18-21 | |
| M | 8.4-9.6 | 21-24 | |
| L | 9.6-10.8 | 24-27 | |
| XL | 10.8-12 | 27-30 | |



AD-402



ANKLET

Provides compression for the foot and ankle joint. Spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use.

INDICATIONS:

- Lymphoedema.
Lipedema.
Chronic venous insufficiency.
Varicose veins.
- Post Thrombotic Syndrome.
Venous Stasis Ulcer.
Post Sclerotherapy.
Lipodermatosclerosis.

| SIZE | INCHES | CM | |
|------|----------|-------|--|
| S | 7.2-8.4 | 18-21 | |
| M | 8.4-9.6 | 21-24 | |
| L | 9.6-10.8 | 24-27 | |
| XL | 10.8-12 | 27-30 | |



AD-403



KNEE CAP

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use.

INDICATIONS:

- Arthritic pain.
Osteoarthritis of the knee.
Ligaments or tendon injury.
- Quad cramps.
Knee sprain.
Weak knees

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |

AD-404



KNEE CAP (WITH OPEN PATELLA)

Spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use. Anterior patellar opening with reinforced padding relieves patellar pressure.

INDICATIONS:

Arthritic pain.
Osteoarthritis of the knee.
Ligaments or tendon injury.

Quad cramps.
Knee sprain.
Weak knees

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-405



KNEE CAP (WITH RIGID HINGE)

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use. Rigid hinges provide extra support.

INDICATIONS:

Arthritic pain.
Osteoarthritis of the knee.
Ligaments or tendon injury.

Quad cramps.
Knee Sprain.
Weak knees

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-415



KNEE WRAP (NEOPRENE)

Made with triple-layer neoprene fabric. The triple hook loop closing system ensures easy wearing, removal and proper fit. Anterior patellar opening with reinforcement relieves patellar pressure and holds the patella in position.

INDICATIONS:

- Arthritic pain.
- Quad cramps.
- Osteoarthritis of the knee.
- Knee sprain.
- Ligaments or tendon injury.
- Weak knees.

| SIZE | INCHES | CM | |
|-----------|---------|-------|--|
| Universal | 14.8-22 | 37-55 | |
| Special | 22-29.2 | 55-73 | |

AD-410



KNEE WRAP HINGED (NEOPRENE)

Simple and quick put on with an easy wrap-around application made from tapered, four-way stretchable neoprene material. Open knee design mitigates pressure from the patella with dual adjustable straps for compression above and below the knee. Biaxial hinges provide support with even flexion and extension movement of the knee.

INDICATIONS:

- Arthritic pain.
- Quad cramps.
- Osteoarthritis of the knee.
- Knee sprain.
- Ligaments or tendon injury.
- Weak knees

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |

AD-409



ELASTIC KNEE SUPPORT

Made from elastic, comfortable and durable material with adjustable strap for proper compression and customized fit. It winds around the knee to support an appropriate flexion and extension movement. Great for supporting an active lifestyle and sports.

INDICATIONS:

Arthritic pain. Quad cramps.
Osteoarthritis of the knee. Knee Sprain.
Ligaments or tendon injury. Weak knees

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-414



PATELLAR SUPPORT

Provides enhanced kinesthesia for promoting blood circulation and healing with customized compression and padding.

INDICATIONS:

Patellofemoral pain syndrome.
Patellar tendonitis and arthritis.
Osgood-Schlatter disease.
Patellar dislocation and subluxation.
Mild instabilities of the knee joint.

| SIZE | INCHES | CM | |
|-----------|--------|-------|--|
| Universal | 12-20 | 30-50 | |

AD-704



COMPRESSION STOCKING BELOW KNEE

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use.

INDICATIONS:

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 9.8-11.3 | 25-29 | |
| M | 11.3-12.8 | 29-33 | |
| L | 12.8-14.3 | 33-36 | |
| XL | 14.3-15.8 | 36-40 | |
| XXL | 15.8-17.3 | 40-44 | |



AD-703



COMPRESSION STOCKING MID THIGH

Anatomically tapered form and spandex-based fabric ensure effective compression with stretchability and comfort. It is made of skin-friendly material providing proper ventilation and comfort for wearing under clothing. High denier count provides durability and sturdiness for extended use.

INDICATIONS:

- Varicose veins.
- Spider veins.
- Fatigue.
- Inflammation of the veins.
- Aching or swollen legs.

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |

AD-406




KNEE IMMOBILIZER (SHORT)

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

INDICATIONS:

- Ligament injuries.
- Acute sprains/strains of the knee.
- Tibial plateau fractures.
- post-cast care.
- Osteochondral injury.
- Post/Preoperative care.
- Patellar tendon injury.

| SIZE | INCHES | CM |  |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-407




KNEE IMMOBILIZER (LONG)

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

INDICATIONS:

- Ligament injuries.
- Acute sprains/strains of the knee.
- Tibial plateau fractures.
- post-cast care.
- Osteochondral injury.
- Post/Preoperative care.
- Patellar tendon injury.

| SIZE | INCHES | CM |  |
|------|-----------|-------|---|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-408




KNEE IMMOBILIZER (22 INCHES)

Hook & loop straps with loop closures and extra strap at the patellar region for proper fitting along with anatomically shaped posterior and lateral splints provide precise immobilization of the knee.

INDICATIONS:

- Ligament injuries.
- Acute sprains/strains of the knee.
- Tibial plateau fractures.
- post-cast care.
- Osteochondral injury.
- Post/Preoperative care.
- Patellar tendon injury.

| SIZE | INCHES | CM |  |
|------|-----------|-------|--|
| S | 14.8-17.2 | 37-43 | |
| M | 17.2-19.6 | 43-49 | |
| L | 19.6 - 22 | 49-55 | |
| XL | 22-24.4 | 55-61 | |
| XXL | 24.4-26.8 | 61-67 | |



AD-411




R.O.M. KNEE BRACE

Our R.O.M knee brace provides locking at any angle or motion between any two angles with intervals of 15 degrees. With neoprene lining and silicone coating it renders cushioning and comfort. It is easy to wear and remove with straps and reverse buckle mechanism which provide proper fitting.

Dial - 0-120° Flexion, 0-60° Extension

INDICATIONS:

- Arthritis/Osteoarthritis.
- Lateral Collateral Swelling & inflammation.
- Hamstring Injury.
- Total knee replacement.
- Knee ligament instability.
- Knee sprain.

| SIZE | INCHES | CM |  |
|-----------|---------|-------|---|
| Universal | 14.8-22 | 37-55 | |

AD-417



CALF SUPPORT

Hypoallergenic elastic fabric from unique fibres blended with cotton prevents skin irritations and perfect compression. Porous weaving allows airflow and proper ventilation keeping the skin dry.

INDICATIONS:

- Weak or injured calves, fatigued legs, and stressed muscles.
- Aching calf muscles.
- Cramps and painful calf conditions.
- Muscle strain.

| SIZE | INCHES | CM | |
|------|-----------|-------|--|
| S | 9.8-11.3 | 25-29 | |
| M | 11.3-12.8 | 29-33 | |
| L | 12.8-14.3 | 33-36 | |
| XL | 14.3-15.8 | 36-40 | |







FOOT, HEEL & TOES SUPPORT

FOOT DROP SPLINT (LEFT/RIGHT)

ARCH SUPPORT

BUNION SPLINT

TOE SEPARATOR

INSOLE (SILICONE)

HEEL CUSHION (SILICONE)

METATARSAL PAD

CAST SHOE



AD-413




FOOT DROP SPLINT

Useful for patients under critical care of foot injury. Lightweight, robust and flexible polypropylene shell in one piece. Comfortable padded calf and ankle straps for additional comfort and fitting.

INDICATIONS:

- Muscle weakness due to stroke, spinal cord injury, muscular dystrophy, cerebral paralysis, poliomyelitis, multiple sclerosis and peripheral neuropathy.
- Recovery from surgical or neurological treatments for foot drop.
Muscle tightening and runner's cramp.
Plantar fasciitis.
Achilles tendonitis.

| SIZE | EUROPE | USA | MM |  |
|------|--------|-----------|---------|--|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |



AD-706




ARCH SUPPORT

Silicone based material absorbs shock, reduces friction and improves general foot comfort. Controls foot position and prevents arch pain. Can be used with socks.

INDICATIONS:

- Flat foot and associated problems.
Walking discomfort.
Plantar fasciitis.

| SIZE | EUROPE | USA | MM |  |
|------|--------|-----------|---------|---|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |

AD-708




BUNION SPLINT

The anatomically moulded splint is cast from high-quality polymer for rigidity and durability. Fastened metal hinges for improved durability. Polyethylene foam provides good cushioning and added comfort. Washable, ventilated and adapts to nearly any size to easily wear with shoes.

INDICATIONS:

Hallux Valgus.
Toe discomfort.

| SIZE |  |
|-----------|--|
| Universal | |



AD-705




TOE SEPARATOR

Odour-free, silicone based toe separator helps release pressure and realign the toes.

INDICATIONS:

Bunion, corns and calluses pain.
Hallux Valgus.
Toe misalignment.
Plantar Fasciitis.
Crooked toes, claw toes & hammertoes.

| SIZE | EUROPE | USA | MM |  |
|------|--------|-----------|---------|---|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |



AD-702



INSOLE

Made of medical grade silicone for durability and comfort. The silicone guarantees an optimal absorption of the shock generated by the daily activities and the high impacts during sports. The softer blue insert relieves the strain on sore heel areas and toes.

INDICATIONS:

- Achilles Tendonitis.
- Heel pain.
- Heel spurs.
- Pain due to walking.

| SIZE | EUROPE | USA | MM | |
|------|--------|-----------|---------|--|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |

AD-701



HEEL CUSHION

Made of medical grade silicone for durability and comfort. The silicone guarantees an optimal absorption of the shock generated by the daily activities and the high impacts during sports. The softer blue insert relieves the strain on sore heel areas.

INDICATIONS:

- Achilles Tendonitis.
- Heel pain.
- Heel spurs.

| SIZE | EUROPE | USA | MM | |
|------|--------|-----------|---------|--|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |

AD-707



METATARSAL PAD

Decreasing the pressure on the ball of your foot by supporting the metatarsal bone. Spreads some of the pressure on the shafts of the bones thereby reducing the pressure on the balls of your feet.

INDICATIONS:

- Plantar Fasciitis.
- Plantar Fat Pad Atrophy.
- Metatarsalgia (stone bruise)

| SIZE |
|-----------|
| Universal |



AD-418



CAST SHOE

A lightweight cast shoe which improves balance and helps early ambulation. Offers cast protection and increases ease of walking. Hook and loop fastener system for a good grip.

INDICATIONS:

- Neuropathic diabetic foot ulcers.
- Forefoot and mid-foot injury.
- Post ligament, soft tissue and tendon surgery.

| SIZE | EUROPE | USA | MM | |
|------|--------|-----------|---------|--|
| S | 34-36 | 2.75-4.35 | 217-230 | |
| M | 37-39 | 5.25-6.75 | 237-250 | |
| L | 40-42 | 7.5-9.0 | 250-270 | |
| XL | 43-44 | 9.5-11 | 277-290 | |





HOT & COLD THERAPY AIDS

ADORE
Supporting Health



AD-1003



COOL EYE PACK

Easy to use eye pack free of bacteria and odours. Made with non-toxic materials and designed to be reused.

INDICATIONS:

For weary and watery eyes.
Useful in redness, burning eyes and conjunctivitis.
Relaxes eyes in insomnia.

SIZE: Universal



AD-1001



HOT & COLD PACK

Pack made of heating/refrigerant gel which instantly relieves muscle pain and other aches due to wounds, sprains, etc. Easy to apply and use.

INDICATIONS:

Lumbago (Low Back Pain).
Joint pain & swelling.
Muscle spasms.

SIZE: Universal



Adore Medical Private Limited

ADORE

Supporting Health



PRODUCT CATALOG 2021



ADORE

Supporting Health



HEAD OFFICE

Plot No. A-189(A),
Road No.01-D,
Vishwakarma Industrial Area,
Jaipur - 302013,
Rajasthan, India

CONTACT SALES

Phone : +91141401 5880
sales@adoremedical.com

ONLINE

www.AdoreMedical.com

